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## User Manual **EN** Digital Pocket Scale F / FH

Thank you for purchasing our product. Please read all operating instructions carefully before use.

- Never load the scale with more than the maximal capacity. Overloading will permanently damage it!
- Avoiding any exposure to extreme heat or cold, your scale works better when operated at normal room temperature. Always allow the scale to acclimate to normal room temperature for at least one hour before use.
- Keep your scale in a clean environment. Dust, dirt, moisture, vibration, air currents and a close distance to other electronic equipment can all cause an adverse effect on the reliability and accuracy of your scale.
- Handle with care. Avoid shaking, dropping the scale.
- Use the scale on a hard flat surface.

### I. Technical data

Model	Min- / Maximum Load	Readability
F-1000	0.5g / 1000g	0.05g
F-2000	1g / 2000g	0.1g
F-3000	1g / 3000g	0.1g
FH-500	0.1g / 500g	0.01g to 200g, 0.05g to 500g
FH-2000	0.5g / 2000g	0.05g to 1000g, 0.1g to 2000g

Weighing units:	g, oz, dwt, ozt, ct, gn
Keyboard:	MODE, TARE, ON/OFF
Power supply:	4 x 1.5V AAA Batteries
Weighing plate (Ø / W x D):	93 x 68 mm
Dimensions (B x T x H):	115 x 90 x 22 mm
Net weight:	210 g
Operating temperature:	10° to 30°C, unto 80% R.H.
Display:	LCD display with backlight, Digit height 10 mm
Auto power off:	After about 60sec. without weight change

## II. Operation

- Place the scale on a horizontal flat surface. Press [ON/OFF] key.
- Wait until the LCD displays "0.00".
- Put the object(s) on the weighing platform.
- Using the [MODE]-key to switch the weighing units.

## III. Tare Weighing

- Turn on scale as described above.
- Place the "tare item" on the platform.
- Press the [TARE] key and wait until "0.00" is displayed.
- Add the "net-weight-item".

## IV. Weighing units

- Turn on scale as described above.
- Press the [MODE] repeatedly to change the weight unit until the required weight unit appears

## V. Calibration

**NOTE:** The scale has been calibrated by the factory. For the calibration, a preset adjustment weight is required. Without such an adjustment weight, calibration is not possible and leads to erroneous values.

### **Please ask the merchant for the suitable test weights.**

- In the switch off state, press the [MODE] key and hold it and then press the [ON/OFF] key without a weight on the scale. Wait until the display shows "888888".
- Then release the buttons. After a few seconds flashes "EXXX" on the display, the Model of the scale: about E1000 for the FH-1000. "E" means the calibrations function is activated.
- Place **the suitable weight** on the weighing surface and wait until the scale indicates "PASS" or "good".
- The process is complete and the value was stored.

## VI. Error Messages and Troubleshooting

"E" or no display	problem of energy supply, check the voltage of batteries and replace them if necessary.
"H"	Capacity exceeded, remove the weights immediately to prevent permanent damage on the scales.
Fluctuating readings	check the environmental conditions.
Incorrect readings	recalibrate the scale is required.

### **G&G GmbH**

Novesiastr. 31  
41564 Kaarst  
Germany

Tel: +49 (0 21 31) 1 33 56 05

Fax: +49 (0 21 31) 1 33 56 07

Email: [info@gandg.de](mailto:info@gandg.de)

Homepage: [www.gandg.de](http://www.gandg.de)