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Digitale Feinwaage TS

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User Manual **EN**

Digital Pocket Scale TS

Thank you for purchasing our product. Please read all operating instructions carefully before use.

- Never load the scale with more than the maximal capacity. Overloading will permanently damage it!
- Avoiding any exposure to extreme heat or cold, your scale works better when operated at normal room temperature. Always allow the scale to acclimate to normal room temperature for at least one hour before use.
- Keep your scale in a clean environment. Dust, dirt, moisture, vibration, air currents and a close distance to other electronic equipment can all cause an adverse effect on the reliability and accuracy of your scale.
- Handle with care. Avoid shaking, dropping the scale.
- Use the scale on a hard flat surface.

I. Technical data

Model	Min- / Maximum Load	Readability
TS-100	0.2g / 100g	0.01g
TS-200	0.2g / 200g	0.01g
TS-300	0.2g / 300g	0.01g
TS-500	0.2g / 500g	0.01g
TS-1000	1g / 1000g	0.1g

Weighing units:	g, oz, dwt, ozt, ct
Keyboard:	ON/OFF, Z/T, MODE, C
Power supply:	2 x AAA Batteries
Weighing plate (Ø / W x D):	60 x 60 mm
Dimensions (B x T x H):	115 x 75 x 20 mm
Net weight:	110 g
Operating temperature:	10° to 30°C, 80% R.H.
Display:	Touch Screen, LCD display with backlight, digit height 12 mm

Auto power off:

After about 60sec. without weight change

II. Operation

- Place the scale on a horizontal flat surface. Press [ON/OFF] key.
- Wait until the LCD displays "0.00".
- Put the object(s) on the weighing platform.
- Using the [MODE]-key to switch the weighing units.

III. Tare Weighing

- Turn on scale as described above.
- Place the "tare item" on the platform.
- Press the [Z/T] key and wait until "0.00" is displayed.
- Add the "net-weight-item".

IV. Weighing units

- Turn on scale as described above.
- Press the [MODE] repeatedly to change the weight unit until the required weight unit appears

V. Counting Function

- Press and hold the [C] key till display shows "PCS". Press the [MODE] key to select a reference amount (25, 50, 75 or 100 pcs)
- Place a certain amount (25, 50, 75 or 100 pcs) of the item which you want to weigh and count onto the platform. Press the [C] key to active the counting function. The LCD will display "XX", and the reference value be confirmed. The counting function is ready!
- With parts counting you can either count parts into a container or remove parts from a container. *As a rule:* The higher the reference quantity the higher the counting exactness
- Now with pressing the [MODE] key to return the weight mode, and press the [C] key to enter the counting mode.

VI. Calibration

NOTE: The scale has been calibrated by the factory. For the calibration, a preset adjustment weight is required. Without such a adjustment weight the calibration is not possible and leads to erroneous values.

Please ask the merchant for the suitable test weights.

- Turn on the scale by pressing the [ON/OFF] without a weight on the scale. Wait until the display shows "0.00".
- Press and hold the [MODE] key for 3 seconds until the display shows "CAL", then release the button. Depending on the model the Display will show "CAL" flashing and automatically continue at the next step, if not press [MODE] again. The Display will now start to flash.
- After a few seconds the Display flashes "XXX g" on the display (required adjustment weight).
- Place the requested weight on the surface and wait until the scale indicates the value. Depending on the model the scale now will show "PASS" or a new calibration-value. If a new weight is required remove the old weight and place the new weight on the surface, wait until it shows "PASS"
- As soon as you see "PASS" the calibration is completed successfully.

G&G GmbH

Novesiastr. 31

41564 Kaarst

Germany

Tel: +49 (0 21 31) 1 33 56 05

Fax: +49 (0 21 31) 1 33 56 07

Email: info@gandg.de

Homepage: www.gandg.de