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Digitale Feinwaage MS

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User Manual **EN**

Digital Pocket Scale MS

Thank you for purchasing our product. Please read all operating instructions carefully before use.

- Never load the scale with more than the maximal capacity. Overloading will permanently damage it!
- Avoiding any exposure to extreme heat or cold, your scale works better when operated at normal room temperature. Always allow the scale to acclimate to normal room temperature for at least one hour before use.
- Keep your scale in a clean environment. Dust, dirt, moisture, vibration, air currents and a close distance to other electronic equipment can all cause an adverse effect on the reliability and accuracy of your scale.
- Handle with care. Avoid shaking, dropping the scale.
- Use the scale on a hard flat surface.

I. Technical data

Model	Min- / Maximum Load	Readability
MS-50	0,1g / 50g	0,005g
MS-100	0,2g / 100g	0,01g
MS-200	0,2g / 200g	0,01g
MS -300	0,2g / 300g	0,01g
MS -500	0,2g / 500g	0,01g
MS -1000	1g / 1000g	0,1g

Weighing units:	g, oz, dwt, ozt, ct
Keyboard:	ON/OFF, TARE, MODE
Power supply:	2 x AAA Batteries
Weighing plate (Ø / W x D):	75 x 63 mm
Dimensions (B x T x H):	128 x 78 x 19 mm
Net weight:	145 g
Operating temperature:	10° to 30°C, 80% R.H.

Display:	LCD display with backlight, digit height 10 mm
Auto power off:	After about 60sec. without weight change

II. Operation

- Place the scale on a horizontal flat surface. Press [ON/OFF] key.
- Wait until the LCD displays "0.00".
- Put the object(s) on the weighing platform.

III. Tare Weighing

- Turn on scale as described above.
- Place the "tare item" on the platform.
- Press the [TARE] key and wait until "0.00" is displayed.
- Add the "net-weight-item".

IV. Weighing units

- Turn on scale as described above.
- Press the [MODE] repeatedly to change the weight unit until the required weight unit appears

V. Calibration

NOTE: The scale has been calibrated by the factory. For the calibration, a preset adjustment weight is required. Without such a adjustment weight the calibration is not possible and leads to erroneous values.

Please ask the merchant for the suitable test weights.

- Turn on the scale by pressing the [ON/OFF] without a weight on the scale. Wait until the display shows "0.00".
- Press and hold the [MODE] key for 3 seconds until the display shows "10000" (or a similar number), then release the button and press it again. The Display will now show "CAL".
- After a few seconds the Display flashes "XXX g" on the display (required adjustment weight).

- Place the requested weight on the surface and wait until the scale indicates the value. Depending on the model the scale now will show “PASS” or a new calibration-value. If a new weight is required remove the old weight and place the new weight on the surface, wait until it shows “PASS”
- As soon as you see “PASS” the calibration is completed successfully.

G&G GmbH

Novesiastr. 31

41564 Kaarst

Germany

Tel: +49 (0 21 31) 1 33 56 05

Fax: +49 (0 21 31) 1 33 56 07

Email: info@gandg.de

Homepage: www.gandg.de